

Protocol for Ussing Chamber Recording

1. Turn on 1) the water bath (42 °C) (NESLAB EX7, Thermo Fisher Scientific, Newington, NH, USA), 2) multichannel voltage/current clamp VCC MC8 (Physiologic Instrument of Medical Research, San Diego, CA, USA), and 3) computer.
2. Preheat all test solutions to 37°C in a water bath.
3. Connect electrodes (two black and two white for each chamber, filled with 4% agarose/3M KCl in the tip and backfilled with 3M KCl) to each chamber in left-to-right order white, black, black, white.
4. Add 5ml bath solution to each side of a chamber (get rid of air bubbles !!!).
5. Turn on the gas cylinder (mixed gas of 5% CO₂ and 95% O₂) and connect the air supply to each chamber half.
6. CORRECT ELECTRODE OFFSET POTENTIAL: Press the MASTER OVERRIDE switch to "ON", the FUNCTION key to "OPEN" and METER switch to "V" to test all 8 channels. . The panel meter will display the difference potential between the pairs of voltage-sensing (black) electrodes. Press the OFFSET button to to the polarity opposite that displayed on the meter and use the dial to adjust the offset so that the meter reads 0.0.
7. LIQUID RESISTANCE CORRECTION: Switch the meter to I. Press the PUSH TO ADJ button to pass current across the fluid in the chamber. The meter should read ~64µA. Switch the METER to V. Press and hold the PUSH TO ADJ button and adjust the FLUID RES COMP so that the displayed voltage reads zero. The resistance of the fluid in the chamber can be read from the fluid resistance dial. This value depends on the slider used and temperature of the bathing solution. For the P2306 slider at 37°C the resistance should be about 540 Ohm).
8. MOUNTING TISSUE PATCH/MONOLAYERS:
Working with one chamber at a time, perform the following:
 - a. Recheck offset and fluid resistance just prior to removing the slider to mount the tissue. Readjust if necessary.
 - b. Remove the air tubes from the chambers;
 - c. Suck the solutions from each side, loosen thumbscrew and spread chambers apart.
 - d. Remove the slider, dry and replace layer of silicone vacuum grease. Remove any saline in the space between the chamber halves.
 - e. Place the tissue over the aperture in the slider;
 - f. Insert the slider containing the tissue or cell culture into the chamber so that the apical side of the cell faces to the left. Ensure that no saline reaches the blue heatblock as this can cause crosstalk between chambers.

- a. Wash the chamber with dH₂O after each experiment. Suck off the water with tissue paper.
 - b. Immerse the tips of the electrodes in physiological saline (sans glucose) or 3M KCl. (I think the saline will provide for less electrode drift because the diffusion of KCl into the saline will be in a steady-state).
19. Turn off power to the amplifier and water bath. Close the main valve on the gas cylinder.
 20. Transfer data to Excel or OriginLab for analysis